# Air Cooling System

# The most powerful Air Cooling System



Model : Cryo - Q

#### ★Convenient and economical for users

- It is very economical because treatment times can be considerably reduced.
- · It is very convenient because using is very easy.
- $\cdot\,$  It is very useful because consumable parts are not needed.

## ★ Very effective and safe for patients

- · It can reduce skin damage because it can be used with lasers.
- · It can reduce pain and recovery time.
- · It can protect skin because cooling can be made evenly.

## ★The suitable lasers

- · Fractional CO2 Laser, IPL, RF Microneedle type, Er YAG,
- · Diode laser, Alexandrite, Nd:YAG(Q-Switched, long pulsed yag)

#### ★ Special Features.

- · Air flow controlling system
- · Cooling air temperature (Max. 35 degree)
- · Quick cooling with high quality compressor
- · Auto adjustable controlling system
- · Auto detecting system to keep best condition



Technical Specifications	
A.C Input	220V,50/60Hz
Air level	10step
Size	400(W)x540(L)x940(H)
Air flow output	1000L/Min
Temperature	0~ - 35°C
Weight	50kg





A.M.I Inc.