



Detail

| | | |
|---------------|--------------------------|---------------|
| Input | AC 100V/240V 50Hz 5~8.5A | |
| Output | DC 48V 13A | |
| Control Board | MAIN BOARD | |
| | LED CONTROL BOARD | |
| | 4 COLOR LED BOARD | |
| Wave | 415nm | Visible rays |
| | 585nm | Visible rays |
| | 635nm | Visible rays |
| | 830nm | Infrared rays |

DR.LUX

The advanced LED light therapy

01 Maximum Effectiveness through Combination of Various Wavelengths


- **Skin Brightening Treatment**
Blue Light (415nm) + Amber Light (585nm) + Red Light (635nm)
 - **Acne care**
Blue Light (415nm) + Red Light (635nm)
 - **Pain Relief**
Red Light (635nm) + Infrared Light (830nm)
 - **Scalp Treatment**
Red Light (635nm) + Infrared Light (830nm)
 - **Redness**
Amber Light (585nm) + Red Light (635nm)
- * Up to 3 wavelengths can be combined.



02 No LED Panel Change needed for Use of Different Light Mode

It is easy to operate in one LED Panel and incorporate into any practices. Four treatments and multiple wavelengths in one panel platform stimulate skin cells in various layers to counteract the effects of aging and environment.

03 Treatment



The advanced LED light therapy technology starts from Dr.Lux. Now, LED gets more values over light with Dr.Lux.

04 4 Wavelengths in 1 System

Promotes calming and relaxation of skin
Assists in killing bacteria in the skin
Helps to normalize oil production
Has a deep cleansing effect
Shrinks enlarged pores
Reduces active acne
Prevents breakouts



Stimulates the nervous and lymphatic system
Decreases skin redness & flushing
Regenerates damaged tissues
Decreases skin irritation
Decreases rosacea



Stimulates the production of collagen and elastin
Reduces the appearance of fine lines, wrinkles and superficial hyperpigmentation
Allows for bald head and psoriasis treatment
Increases circulation and lymph drainage
Improves skin elasticity
Smooths skin texture



Provides anti-aging benefits to your body
Stimulates fibroblasts proliferation
Promotes cell regeneration
Improves blood circulation
Rejuvenates the skin
Heals wounds faster
Provides pain relief

